

WIRELESS AUTO-TERSUS



- Wireless Auto-Tersus is Status Medical Equipment's state of the art Uroflowmeter.
- Undoubtedly the most Hassle free Uroflowmeter in the market.

Advantages over Conventional Uroflowmeter

Autoflush cleans the Urinal after every test is completed, No need for the clinical staff to clean the beaker after test.

Male Patient voids in the most natural position avoiding all the artifacts resulting from awkwardness and occasional spilling in the conventional uroflowmeter.

No need for the wired connection from the bathroom to the computer.

Wireless connection allows the urinal to be in a bathroom whereas all the reports can be monitored and printed from computer in attendant's cabin.

Attendant can Add patient, Start the test, Stop the test, Print the report. All this while sitting in his cabin using a StatFlow®, Status Medical Equipment's proprietary software.

System Consists of:


- Auto-Tersus Urinal
- Control Unit
- Processing Unit with pre- installed StatFlow® OS
- Transreceiver module




Metro Health Hospital

4th Avenue, Centre Street
New York, USA

Patient Name: Andrew K Ruffalo
Patient ID: MHH112 Test Date: 02/Aug/2017-15:47 Gender: Male
Height: 150 cm Weight: 52 kg Age: 17 years
Position: Standing Referred By: Ref Attended By: Att



| | | | | | |
|-------------------|-----------|-----------------|-----------|---------------------------|--------|
| Hesitancy: | 9 sec | Flow time: | 24 sec | Voided Volume: | 330 ml |
| Time to max flow: | 12 sec | Max Flow: | 20 ml/sec | Voiding time: | 24 sec |
| Avg Flow: | 14 ml/sec | Residual urine: | - ml | Post test residual urine: | 30 ml |



ASTRA 30


EVEREST

URODYNAMICS

UROFLOWMETER

WORK STATION

HYDROMAT



STATUS MEDICAL EQUIPMENTS

M-9/A, Additional M.I.D.C. Satara 415004,
Maharashtra India

Phone 91-02162-240037, 240045,

Fax 91-02162-240743

e-mail statmedi@gmail.com

www.statmedi.com