

# Wireless FlowComp-NG



## The wireless uroflowmeter



- Status Medical Equipment presents Wireless FlowComp-NG, the wireless uroflowmeter.
- The RF based wireless technology of FlowComp-NG enables the fast and flawless communication.

### Advantages over the Conventional Uroflowmeter:

- No need for the wired connection from bathroom to computer.
- Wireless connection allows the Vmic sensor to be kept in the Bathroom whereas the reports can be monitored from computer in attendant's cabin.
- FlowComp-NG comes with the StatFlow, SME's proprietary software specially developed for uroflowmetry.
- Using the StatFlow test attendant can add patient, start test, stop test, print the report while sitting in his cabin.
- You will not need a separate computer for uroflowmetry, only the LCD screen is sufficient.
- Sampling rate, Reporting parameters and accuracy of the readings is as per the ICS (International Continence Society) guidelines.
- The new report format is well organised and easy to read consisting all the report parameters required for the uroflowmetry diagnosis.
- As the FlowComp-NG is computer based, SME can provide the online support for this machine which will help us to serve you better and quick.



### Metro Health Hospital

4th Avenue, Centre Street  
New York, USA

Patient Name: Andrew K Ruffalo

Patient ID: MHH112 Test Date: 02/Aug/2017-15:47


Height: 150 cm Weight: 52 kg

Position: Standing Referred By: Ref


Gender: Male

Age: 17 years

Attended By: Att



|                   |           |                 |           |                           |        |
|-------------------|-----------|-----------------|-----------|---------------------------|--------|
| Hesitancy:        | 9 sec     | Flow time:      | 24 sec    | Voided Volume:            | 330 ml |
| Time to max flow: | 12 sec    | Max Flow:       | 20 ml/sec | Voiding time:             | 24 sec |
| Avg Flow:         | 14 ml/sec | Residual urine: | - ml      | Post test residual urine: | 30 ml  |



ASTRA 30


EVEREST

URODYNAMICS

UROFLOWMETER

WORK STATION

HYDROMAT



## STATUS MEDICAL EQUIPMENTS

M-9/A, Additional M.I.D.C. Satara 415004,  
Maharashtra India

Phone 91-02162-240037, 240045,

Fax 91-02162-240743

e-mail [statmedi@gmail.com](mailto:statmedi@gmail.com)

[www.statmedi.com](http://www.statmedi.com)